

Micromotor for Laboratory Work

FORZA L50K

OPERATION MANUAL




OM-E387E

Please read this operation manual carefully before use and keep for future reference.

Thank you very much for purchasing the Forza L50K.

Please read this Operation Manual carefully regarding the instructions for use, handling method, and maintenance and inspection prior to use. Store the manual in a place where users can review it anytime.

Contents

 IMPORTANT INSTRUCTIONS AND WARNING	1
1. Component Names	5
2. Set up of Control Unit	6
3. Operation Procedure	7
4. Speed limit mechanism.....	8
5. Protective Circuit for Motor	8
6. Memory Function	8
7. Error Code	9
8. Replacement of Fuse	10
9. Maintenance Mode	10
10. Vacuum-coupled Mode	11
11. Handling of Motor and Handpiece	12
12. Handpiece Stand	14
13. Handpiece Holder	14
14. Specification	14
15. Troubleshooting	15
16. Disposing Product	17
17. Warranty	17

IMPORTANT INSTRUCTIONS AND WARNING -Electric Devices




WARNING!

When using electric tools, basic safety precautions should always be followed to reduce the risk of fire, electrical shock and personal injury. Read all these instructions before operating this product and save these instructions.

A. GROUNDING INSTRUCTIONS

1. In the event of a malfunction or breakdown, grounding provides a path of least resistance for electric current to reduce the risk of electric shock. This tool is equipped with an electric cord having an equipment-grounding conductor and a grounding plug. The plug must be plugged into a matching outlet that is properly installed and grounded in accordance with all local codes and ordinances.
2. Do not modify the plug provided - if it will not fit in the outlet, have the proper outlet installed by a qualified electrician.
3. Improper connection of the equipment-grounding conductor can result in a risk of electric shock. The conductor with insulation having an outer surface that is green with or without yellow stripes is the equipment-grounding conductor. If repair or replacement of the electric cord or plug is necessary, do not connect the equipment grounding conductor to a live terminal.
4. Check with a qualified electrician or service personnel if the grounding instructions are not completely understood, or if in doubt as to whether the tool is properly grounded.
5. Use only 3-wire extension cords that have 3-prong grounding plugs and 3-pole receptacles that accept the tool's plug.
6. Repair or replace damaged or worn cord immediately.
7. This tool is intended for use on a circuit that has an outlet that looks like the one illustrated in Sketch A in Figure (below) (120V). The tool has a grounding plug that looks like the plug illustrated in Sketch A in Figure (below). A temporary adapter, which looks like the adapter illustrated in Sketches B and C, may be used to connect this plug to a 2-pole receptacle as shown in Sketch B if a properly grounded outlet is not available. The temporary adapter should be used only until a properly grounded outlet can be installed by a qualified electrician. The green-colored rigid ear, lug, and the like, extending from the adapter must be connected to a permanent ground such as a properly grounded outlet box.

Grounding Method



8. **USE PROPER EXTENSION CORD.** Make sure your extension cord is in good condition. When using an extension cord, be sure to use one heavy enough to carry the current your product will draw. An undersized cord will cause a drop in line voltage resulting in loss of power and overheating. Table (below) shows the correct size to use depending on cord length and nameplate ampere rating. If in doubt, use the next heavier gauge. The smaller the gauge number, the heavier the cord.
9. Install equipment so that the power supply cord can be pulled out without hindrance in event of emergency. Do not place anything within 10 centimeters around the unit.

Minimum gauge for cord

Ampere Rating		Volts	Total Length of cord			
		120V	7.5m (25ft.)	15m (50ft.)	30m (100ft.)	45m (150ft.)
		240V	15m (50ft.)	30m (100ft.)	60m (200ft.)	90m (300ft.)
More Than	Not More	Code Number				
0	Than		#18	#16	#16	#14
6	6		#18	#16	#14	#12
10	10		#16	#16	#14	#12
12	12		#14	#12	Not Recommended	

B. OTHER WARNING INSTRUCTIONS




1. For your own safety read instruction manual before operating tool.
2. Wear eye protection.
3. Replace cracked disc or polisher immediately.
4. Always use guards and eye shields.
5. Do not over tighten wheel nut.
6. Only use wrench furnished with the Forza L50K.
7. **REMOVE ADJUSTING KEYS AND WRENCHES.** Form a habit of checking to see that keys and adjusting wrenches are removed from tool before turning it on.
8. **KEEP WORK AREA CLEAN.** Cluttered areas and benches invite accidents.
9. **DON'T USE IN DANGEROUS ENVIRONMENT.** Don't use power tools in damp or wet locations, or expose them to rain. Keep work area well lighted.
10. Risk of injury due accidental starting. Do not use in an area where children may be present.
11. **DON'T FORCE TOOL.** It will do the job better and safer at the rate for which it was designed.
12. **USE RIGHT TOOL.** Don't force tool or attachment to do a job for which it was not designed.
13. **WEAR PROPER APPAREL.** Do not wear loose clothing, gloves, neckties, rings, bracelets, or other jewelry that might get caught in moving parts. Nonslip footwear is recommended. Wear protective hair covering to contain long hair.
14. **ALWAYS USE SAFETY GLASSES.** Everyday eyeglasses only have impact resistant lenses, they are NOT safety glasses. Also use face or dust mask if cutting operation is dusty.
15. **SECURE WORK.** Use clamps or a vise to hold work when practical. It's safer than using your hand and it frees both hands to operate tool.
16. **MAINTAIN TOOLS WITH CARE.** Keep tools sharp and clean for best performance and to reduce the risk of injury to persons. Follow instructions for lubricating and changing accessories.

17. DISCONNECT TOOLS before servicing; when changing accessories, such as blades, bits, cutters, and like.
18. REDUCE THE RISK OF UNINTENTIONAL STARTING. Make sure switch is in off position before plugging in.
19. USE RECOMMENDED ACCESSORIES. Consult the owner's manual for recommended accessories. The use of improper accessories may cause risk of injury to persons.
20. NEVER LEAVE TOOL RUNNING UNATTENDED. TURN POWER OFF. Don't leave tool until it comes to a complete stop.
21. For recommended operating speed for various applications, please follow the instructions of bur manufacturers.
22. The system functions normally in the environment where the temperature is at 0-40°C (32-104°F), humidity at 10-85% RH, atmospheric pressure at 500-1060hPa, and no moisture condensation in the Unit. Use at outside of these limits may cause malfunction.
23. Store the system in the place where the temperature is at -10-60°C (14-140°F), humidity at 10-85% RH, atmospheric pressure at 500-1060hPa, and the system is not subject to air with dust, sulfur, or salinity.
24. Severe shock – Eg. Dropping Control Unit, or the Micromotor – may cause damage.
25. Do not turn the bur lock ring while the handpiece is rotating.
26. Do not rotate the motor when the bur lock ring is at OPEN position, or if a bur is not mounted in the chuck. It may cause the motor/handpiece disconnection or sudden heat generation.
27. Activation of Circuit Breaker means too much load is applied to the motor beyond the motor's capacity. This circuit breaker is designed to protect the motor, but it is desired to perform the grinding work without activating the circuit breaker.
28. Care should be taken not to drop micromotor handpiece on floor or hard work surface in order to avoid damage caused by impact shock.

C. Important Instructions and Warning on Forza L50K.

No lubrication is required to either motor or handpiece because ball bearings are impregnated with grease in both motor and handpiece.

1. Component Names



- ① Control Unit
- ② Motor Handpiece
- ③ Foot Pedal
- ④ Power Switch
- ⑤ Motor Connector
- ⑥ Handpiece Stand
- ⑦ Forward/Reverse Selector Switch
- ⑧ Hand/Foot Selector Switch
- ⑨ Motor Switch
- ⑩ Motor LED
- ⑪ Speed Indicator
- ⑫ Speed Control Knob
- ⑬ Speed Limit Release Button
- ⑭ Foot Pedal Connector
- ⑮ Fuse Box
- ⑯ Power supply cord

2. Set up of Control Unit

1. Mounting of Motor

Insert the motor cord plug locator into ⑤ Motor Connector aligning it with the groove of the connector, and tighten the motor cord by rotating plug nut clockwise. (Fig. 1)



Fig.1

2. Mounting of Foot Pedal

Insert the Foot Pedal cord plug locator into ⑭ Foot Pedal Connector aligning it with the groove of the connector. (Fig. 2)



Fig.2

3. Connecting of Power Cord

Securely insert the plug of ⑯ Power Cord into connector at the back of the unit aligning it with the configuration. (Fig. 3)



Fig.3

3. Operation Procedure

1. ⑩ Connect Power Cord to a wall outlet.
2. Make sure that ⑫ Speed Control Knob is at the lowest speed position.
3. Turn ④ Power Switch ON (Green Power Lamp will light). The preset rotation speed will appear on the Indicator.
4. Select the rotation direction with ⑦ Forward/Reverse Selector Switch. Each time this switch is pressed, the direction changes between FORWARD and REVERSE.
5. Select the control method with ⑧ Hand/Foot Selector Switch. Each time this switch is pressed, the method changes between the HAND and FOOT.

Operation - 1

◆HAND Operation

6. Select HAND with ⑧ Hand/Foot Selector Switch.
7. Preset the rotation speed with ⑫ Speed Control Knob with checking the speed on the Speed Display.
8. Press ⑨ Motor Switch, and the motor will run.
9. To stop the motor, press ⑨ Motor Switch again.

Operation - 2

◆Operation by Foot Pedal

6. Select FOOT with ⑧ Hand/Foot Selector Switch.
7. Preset the maximum rotation speed with ⑫ Speed Control Knob with checking the speed on the Speed Display.
8. Depress ③ Foot Pedal, and the motor will run. The rotation speed can be variably controlled within the preset maximum rotation speed range according to the amount of depressing ③ Foot Pedal.

※ Auto speed mechanism

To Fix the speed within the rotation speed range set by the volume, press ⑨ Motor Switch while the motor is running at the desired speed. The display lamp will flash and the rotation speed can be maintained even if ③ Foot Pedal is released. To cancel it, press ⑨ Motor Switch again or depress ③ Foot Pedal again.

4. Speed limit mechanism

For the use of a small-diameter round bur or fissure bur at $40,000\text{min}^{-1}(\text{rpm})$ or more, pressing the ⑬ Speed Limit Release Button, which is equipped with the ⑫ Speed Control Knob, and turning the ⑫ Speed Control Knob allows a setting of up to $50,000\text{min}^{-1}(\text{rpm})$. (Fig.4)



Fig.4



CAUTION

Use of $40,000\text{min}^{-1}(\text{rpm})$ or more is allowable only when the bur manufacturer specifies that the bur's acceptable rotation speed is $40,000\text{min}^{-1}(\text{rpm})$ or more. If the use of more than the acceptable rotation speed is made, the bur may be broken.

5. Protective Circuit for Motor

When the motor is operated with a load exceeding the limit or the handpiece is in an unrotational condition, the circuit to protect the motor and unit works to stop the power supply to the motor, an error code appears on

⑪ Indicator. (Fig.5)



Fig.5

How to reset the protective circuit

During manual operation, it can be reset by pressing ⑨ Motor Switch again. During operation by Foot Pedal, depress ③ Foot Pedal back to stop, and the protective circuit will be reset.

6. Memory Function

When the Power Switch is turned on, the rotation direction and HAND/FOOT selections made when the main switch was last turned off are restored. Special attention should be given to the rotation direction.

7. Error Code

When the motor is stopped due to some trouble such as a failure, overload, wire breakage or misuse, the Speed Display indicates the error code for checking the status of the unit and understanding the cause of the trouble easily.

Error code	Description	Cause
E0	Self-check error	<ul style="list-style-type: none"> • Abnormal internal memory • Broken internal memory
E1	Overcurrent detection error (Hard)	<ul style="list-style-type: none"> • Long-time use at a high load (overcorrect) • Shorted cord (power line) • Shorted motor winding
E2	Overvoltage detection error	<ul style="list-style-type: none"> • Shorted cord (power line), damaged circuit • Broken internal circuit
E3	Motor sensor error	<ul style="list-style-type: none"> • Faulty sensor (Hall IC) in the motor • Disconnected motor cord • Severed cord (signal line) • Open chuck • Faulty handpiece
E4	Unit overheat error	<ul style="list-style-type: none"> • Temperature rise in the unit due to longtime use at a high load • Unit placed under high temperature
E5	PAM circuit error	<ul style="list-style-type: none"> • Abnormal voltage generated in start / stop circuit • Faulty start / stop circuit from PAM (L Slide)
E6	Rotor lock error	<ul style="list-style-type: none"> • Open chuck • Faulty handpiece • Motor Faulty • Faulty sensor (Hall IC) in the motor • Severed cord (signal, power line)
E8	Overvoltage detection error (Soft)	<ul style="list-style-type: none"> • Long-time use at a high load (overcorrect) • Shorted cord (power line) • Shortstop of the motor winding
E9	ITRIP error	<ul style="list-style-type: none"> • Faulty motor and circuit
EF	Foot pedal error	<ul style="list-style-type: none"> • Faulty Foot Pedal or Shorted Foot Pedal cord • Broken internal circuit

For countermeasures against error displays, see the section on 15.Troubleshooting.

8. Replacement of Fuse

Fuse is located in ⑮ Fuse Box. Release ratchet clamp located on the right-and-left top and bottom of the ⑮ Fuse Box and pull it out to change the fuse (T3.15AL250V). (Fig. 6) Insert new fuses (T3.15AL250V) into the ⑮ Fuse Box and then push the ⑮ Fuse Box into the original position until you hear a click. (Fig.6)




Fig.6

- Fuse: Order No. (120V) : 5018427UO



CAUTION

Fuse is burned out when a short circuit occurs or when over-voltage is flowed into the primary current source. If the cause is uncertain, return the product to Brasseler USA for inspection.

9. Maintenance Mode

The unit is provided with a maintenance mode to check the function of the switches, display, Foot Pedal, motor, etc. While pressing ⑧ Hand/Foot Selector Switch and ⑨ Motor Switch at the same time, turn on the Power Switch and keep pressing the button until beeps are made (for about 2 seconds). At this time, turning the volume from the minimum position in order displays "oP", "dP", "HL", "Pd" and "in", allowing the following checks. To release Maintenance Mode, turn Power Switch off and switch on again.

(1) "oP": Switch check (operation check)

Press the switches on the panel, and the right and/or left lamps will light to check to see if the switches operate normally.


(2) "dP": Display check

Press ⑦ Forward/Reverse Selector Switch, and the lamps will light one by one to check to see if they operate normally. To cancel this check, press ⑦ Forward/Reverse Selector Switch again.

(3) "HL": Motor signal check (Hall IC check)

Press ⑦ Forward/Reverse Selector Switch, and ⑪ Indicator will indicate one or two horizontal lines. Turn the motor slowly by hand, and this display will change to one line, two lines, one line, two lines, smoothly from the top to bottom or from the bottom to top. If any one of the three lines does not light, the sensor (Hall IC) in the motor is faulty or the cord is severed, therefore repair is needed. To cancel this check, press ⑦ Forward/Reverse Selector Switch again.

Eg.



(4) "Pd": Foot Pedal check

Press ⑦ Forward/Reverse Selector Switch, and ⑪ Indicator will change. During normal time, Indicator changes in hexadecimals (0~9, A~F) according to the amount of depressing ③ Foot Pedal . Also, depressing the pedal slightly lights Reset Lamp, and depressing it fully extinguishes the lamp. If ⑪ Indicator does not change smoothly or Reset Lamp does not light properly, ③ Foot Pedal may be faulty. To cancel this check, press ⑦ Forward/Reverse Selector Switch again.

(5) "in": Initializing function

Press ⑦ Forward/Reverse Selector Switch, and beeps will be made and rotation direction, HAND/FOOT and other settings will return to the factory set condition.

Rotation direction: FWD (forward)
HAND/FOOT: HAND
Vacuum-coupled mode: OFF

10. Vacuum-coupled Mode

On some dental benches with vacuum dust collector, the motor may be used while being coupled with a dust collector. When such a dust collector * is used, power consumption of Forza L50K can be regulated so that the vacuum-coupled function can work. If you need to connect to a vacuum dust collector, select the mode as follows:

How to select the mode

While pressing ⑦ Forward/Reverse Selector Switch, turn on ④ Power Switch, and the mode can be selected. A long beep indicates vacuum-coupled mode and 2 short beeps indicate non-coupled (energy-saving) mode.

* Each time the switch selection is made, the mode changes between vacuum-coupled mode and non-coupled mode.

11. Handling of Motor and Handpiece

(1) Insertion or Removal of Bur

The chuck is opened by turning the Bur Lock Ring to an open position. The chuck is loosened and the bur can be removed. By turning the ring in the LOCK direction, the chuck is closed and the bur can be mounted. At this time, turn the ring until it clicks. (Fig.7)




Fig.7



CAUTION

- Do not turn the Ring during rotation.
- Do not operate the handpiece without a bur or if the bur lock ring is in the open position. It may be the reason for the failure or overheating of the handpiece.

(2) Cleaning and Replacement of Chuck

(1) Removal of Chuck

To remove the chuck, open the bur lock ring and turn the chuck counterclockwise with the provided chuck wrench. (Fig.8)




Fig.8

* If a large diameter bur is used under high torque, the chuck may rotate in the close direction and the bur may be stuck and cannot be removed. In this case, align the nose's slit and spindle's spanner position (flat part), and apply an L spanner to fix the spindle. Open the ring and turn the chuck counterclockwise with the provided chuck wrench to remove it.

(Fig. 9)




Fig.9

(2) Cleaning of Chuck

Remove and clean the chuck as frequently as possible in the ultrasonic cleaner. Cleaning once a week is recommended.



CAUTION

Neglecting to clean the chuck for a long time is very dangerous because wax, gypsum, etc., accumulate in the chuck and the bur is insecure, causing runout.

(3) Insertion of Bur into the Chuck

Thinly apply oil before insertion. Open the ring, insert the dummy bur or the bur in use into the chuck, and turn the chuck clockwise by hand until it stops. Then, lock the ring, and the chuck will hold the bur securely.

(Fig. 10)




Fig.10



CAUTION

Before using the handpiece, make sure to pull the rotating instrument (bur, etc.) to check that it is securely mounted.

3. Disconnecting and Connecting of Motor Cord to Motor

Remove the cord nut at the rear end of the motor, and the motor cord connector can be pulled out. For connection, align the connector pin and the hole in the motor cord connector, and insert the connector straight until it stops. Then, tighten the cord nut.

* When inserting the connector, do not turn or twist it.



Fig.11

4. Disconnecting Handpiece from Motor

The handpiece and motor are screwed at the midpoint. Firmly grasp the motor outer case and the handpiece outer case and turn it counterclockwise to disconnect. (Fig. 12)



Fig.12



CAUTION for Handpiece connection

When connecting the handpiece to the motor, turn the handpiece clockwise and tighten firmly. If the clutch is not engaged properly, the handpiece cannot be tightened completely. In such case, Do Not Force.

Loosen the handpiece and turn the bur briefly to re-position the clutch. Reconnect the handpiece and tighten securely. (Fig. 13).




Fig.13

12. Handpiece Stand

At the bottom of the handpiece stand, tools necessary for handpiece maintenance and a spare chuck (optional) can be mounted.

(Fig. 14)



Fig.14

13. Handpiece Holder

The handpiece holder that can be mounted on the right side of the control unit will widen your effective work area. Insert the holder into the slot and fasten the screw. Adjust the holder angle to your best position. (Fig. 15)




Fig.15

14. Specification

Control Unit

Model	NE213-G
Power supply	AC120V 50/60 Hz 46VA
Rated Power	46VA
Weight	2.8kg
Dimensions	W167 x D247 x H77mm

Food Pedal

Model	FC-64
Weight	460g

Handpiece

Model	Handpiece Stand
Weight	120g

Motor Handpiece

Model	UMXL-T(P) / UMXL-T(P) 2M
Speed	1,000~50,000 min ⁻¹
Max. Torque	8.7N · cm
Weight	215g (exclude the cord)
Dimensions	L151 x ø29mm
Cord Length	1.2m / 2.0m

15. Troubleshooting

Please check the following points before sending instruments in for repair.

<Unit / Motor>

Trouble		Cause / Check	Remedy
Pilot Lamp does not light		The power plug is disconnected.	Insert the power plug correctly.
		The fuse is blown.	Replace it with a specified fuse. If the reason the fuse has blown is unknown, ask for an inspection.
		Power Switch is faulty.	Ask for repair.
Control Unit display “—“and “Setting Rotation Speed” alternately.		Power on the Control Unit while pressing the Foot Pedal.	Power on the Control Unit without pressing the Foot Pedal. When it does not work definitely, take off the Foot Pedal cord plug from the Control Unit and Power it on. If it does not operate correctly, contact Brasseler USA. After you tried above and the Control Unit could not operate normally, contact Brasseler USA for repair.
The motor and handpiece do not run, or Reset Lamp lights.	Foot Pedal does not work.	The connection of the foot control cord plug is loose.	Connect the foot control cord plug correctly.
		Hand/Foot Selector Switch is set to HAND MODE.	Set Hand/Foot Selector Switch to FOOT MODE.
		Maintenance mode (4) “Pd”: Check to see if the foot control operates normally.	If the foot control does not operate normally, ask for repair of the foot control or replace it with a new one.
	Error code E0 appears.	Turn on the power again.	If the same error code appears, ask for repair of the unit.
	Error code E1 appears.	Turn on the power again.	If it operates normally, the error display is temporarily due to overload, which is not a problem.
		If you have two or more units, replace the motor and the motor cord and check the operation.	If it operates normally after replacing the motor and the motor cord, the motor and/or the motor cord may be shorted. Ask for repair of the motor and/or the motor cord. If the same error code still appears after replacing, ask for repair of the unit.
	Error code E2 appears.	The motor cord is disconnected.	Connect the motor cord correctly.
		If you have two or more units, replace the motor and the motor cord and check the operation.	If it operates normally after replacing the motor and the motor cord, the motor and/or the motor cord may be severed. Ask for repair of the motor and/or the motor cord. If the same error code still appears after replacing, ask for repair of the unit.

Trouble		Cause / Check	Remedy
The motor and handpiece do not run, or Reset Lamp lights.	Error code E3 appears.	The motor cord is disconnected.	Connect the motor cord correctly.
		Maintenance mode (3) "HL": Confirm whether you work by a motor signal check normally.	If any problem is found during a check, the motor cord may be severed or the sensor in the motor may be faulty. Ask for repair.
		The chuck is open.	Lock the chuck. If any problem is found during a check, the motor cord may be severed or the sensor in the motor may be faulty. Ask for repair.
		Check by turning a bur with a finger lightly and rotating it lightly.	If the rotation is abnormal, ask for repair of the motor and handpiece.
	Error code E4 appears.	Vacuum-coupled Mode is on, while the dust collector is not used.	Make the Vacuum-coupled Mode to be off. See Vacuum-coupled Mode section in the manual. (See * Vacuum-coupled Mode)
		After stopping to cool it down for about 10 minutes, check the operation again.	If it operates normally, there is no problem. Check the operating environment, storage location, etc., for high temperature. If the same error code appears frequently, ask for repair of the unit.
	Error code E5 appears.	Turn on the power again, and repeat starting and stopping several times.	If it operates normally, there is no problem. If the same error code appears, ask for repair of the unit.
	Error code E6 appears.	The chuck is open.	Lock the chuck. If any problem is found during a check, the motor cord may be severed or the sensor in the motor may be faulty. Ask for repair.
		Maintenance mode (3) "HL": Confirm whether you work by a motor signal check normally.	If any problem is found during a check, the motor cord may be severed or the sensor in the motor may be faulty. Ask for repair.
		Check by turning a bur with a finger lightly and it lightly.	If the rotation is abnormal, ask for repair of the motor and handpiece.
	Error code E8 appears.	Turn on the power again.	Error indication is a temporary thing caused by the overload, if it works normally. There is no malfunction.
		If you have two or more units, replace the motor cord and check the operation.	If it operates normally after replacing the motor cord, the motor and/or the motor cord may be shorted. Ask for repair of the motor and/or the motor cord. If the same error code still appears after replacing the motor cord, ask for repair of the unit.

Trouble		Cause / Check	Remedy
The motor and handpiece do not run, or Reset Lamp lights.	Error code E9 appears.	If you have two or more units, replace the motor cord and check the operation.	If it operates normally after replacing the motor cord, the motor and/or the motor cord may be shorted. Ask for repair of the motor and/or the motor cord. If the same error code still appears after replacing the motor cord, ask for repair of the unit.
	Error code EF appears.	Maintenance Mode (4) "Pd": Check whether Foot Pedal operates normally by Foot Pedal Mode.	If it does not operate normally, change the Foot Pedal for new or ask for repair.
The rotation speed does not rise.		The maximum rotation speed for operation by foot pedal should be set with the Speed Control Knob.	Set the maximum rotation speed with the Speed Control Knob.

<Handpiece>

Trouble	Cause / Check	Remedy
The handpiece does not run with the chuck tightened.	Entry of foreign matter in the ball bearings or seizure.	Send it to Brasseler USA. Ask for repair.
Heat is generated during rotation.	Entry of foreign matter in the ball bearings, causing wear of the bearings.	Same as the above.
Vibration or noise occurs during rotation.	Same as the above.	Same as the above.
	A bent bur is used.	Replace the bur
Run out of the bur is heavy.	Dust may be stuck in the chuck or spindle.	Clean the inside of the chuck and spindle.
	The chuck is worn.	Replace the chuck.
	The ball bearings are worn.	Send it to Brasseler USA.
The bur comes out	The chuck is loose.	Tighten the chuck securely. (See *Handling of Motor and Handpiece.)

16. Disposing Product

Please consult with Brasseler USA about waste disposal.

17. Warranty

Brasseler USA warrants the handpiece against defects in manufacturing, workmanship and materials. Brasseler USA reserves the right to analyze and determine the cause of any problem. Warranty is voided should the handpiece not be used in accordance with this manual.