

solutions™

Dialite® L20 “Skinny” Thin Disc System



Order Number:
K0229






Figure Number	L20DCHP	L20DMHP	L20DFHP
Head Size ($\frac{1}{10}$ mm)	200	200	200
Head Length (mm)	1.0	1.0	1.0

Only Available Direct from your Brasseler Representative
To order call 800.841.4522 or fax 888.610.1937
In Canada call 800.363.3838 or fax 866.330.4454
Visit our website at BrasselerUSA.com

solutions™

Dialite® L20 “Skinny” Thin Disc System

	Order #	Grit	Usage
	L20DCHP-HP	C	Polisher - Thin disc for interproximal polishing of tight contacts
	L20DMHP-HP	M	Polisher - Thin disc for interproximal polishing of tight contacts
	L20DFHP-HP	F	Polisher - Thin disc for interproximal polishing of tight contacts

TIPS & TECHNIQUES: Works well for incisal cutbacks.

Complimentary Items: Dialite HP Polishers, Dialite Extra-Oral Diamond Paste (5015241U0), Diamond Discs (911H, 911HH, 911HV)

Power Suggestion: Forza L50K or Forza FIS Electric Attachment

Sterilization: Autoclave only

Recommended Speeds: 8,000-10,000 rpms

By Your Side In Dentistry



BRASSLER USA
DENTAL INSTRUMENTATION