



BRASSELER
USA®

HYGIENE PROMOTIONS | 2015



Prices good through December 31, 2015

BY YOUR SIDE

HANDLE TYPES

#2 Handle



- 1/4 in round, solid stainless steel
- Also known as standard handle
- 26 grams

#4 Handle



- 5/16 in round, hollow stainless steel with matte finish
- Ideal for smaller hands, the hollow core provides ideal tactile sensitivity
- 16.8 grams

#6 Handle




- 3/8 in round hollow stainless steel with matte finish
- Extended handle neck assists in correct finger placement
- 16.4 grams

#8 Cushion Grip Handle



- 7/16 in round with silicone grips around a hollow stainless steel core
- Features an integrated code ring system and elliptically shaped grips to allow customized finger placement
- 20 grams




DIAGNOSTIC

- Probes and explorers manufactured from spring-tempered stainless steel wire for maximum durability and enhanced tactile sensitivity.
- Color-coded probes feature clear, laser-etched markings and will not chip or fade.

► MIRROR HANDLES

25% OFF 6 OR MORE



MH6CS

Lightweight ergonomic diamond knurled, #6 cone socket
• Hollow stainless steel

\$18.50 ea.



MH1CS

Standard handle, #2 cone socket
• Solid stainless steel

\$11.50 ea.

► MIRROR FACES

20% OFF 3 OR MORE BOXES OF 12



FSCS4/1

\$4.16 ea.

FSCS4/12

\$49.95/Box

• Front surface, cone socket
• #4 7/8 in (22 mm) diameter
• Sold individually or in box of 12



FSCS5/1

\$4.16 ea.

FSCS5/12

\$49.95/Box

• Front surface, cone socket
• #5 15/16 in (24 mm) diameter
• Sold individually or in box of 12

► COTTON & DRESSING PLIERS

25% OFF 6 OR MORE



100

Non-Locking Dressing Pliers
• 6 in / 150 mm

\$24.50 ea.



200

Locking Dressing Pliers
• 6 in / 150 mm

\$37.50 ea.



400

Non-Locking Perry Dressing Pliers
• 5 in / 125 mm

\$24.50 ea.

DIAGNOSTIC

DOUBLE-END EXPLORERS

25% OFF 6 OR MORE



EX5

Combination of Explorer 23 and 17
• Instrument of choice for caries and calculus detection

#2, #4, #6

\$18.50 ea.



EX0DU11/12

Old Dominion University 11-12
• Instrument of choice for posterior calculus detection
• One instrument is used for exploring entire dentition

#2, #6

\$20.50 ea.

EXPROS

25% OFF 6 OR MORE



EX23/UNC12

Explorer 23 • Probe end color coded UNC12

#6

\$28.50 ea.



EX23/CP12

Explorer 23 • Probe end color coded Marquis

#2, #4, #6

\$28.50 ea.



EX23/POW


Explorer 23 • Probe end also known as Michigan "O" with Williams markings (Fine Round)

#2, #4, #6

\$28.50 ea.

DOUBLE-END PROBES

25% OFF 6 OR MORE




CP12D

Offset Marquis
• Also known as CP12

#6

\$28.50 ea.




CP12/COW

Marquis/Williams
• Color coded

#6

\$28.50 ea.



CPOWD

Williams Offset
• Color coded

#2, #6


\$28.50 ea.

SCALERS

- Sickle scalers are designed almost exclusively for supragingival periodontal debridement procedures.
- Manufactured from 440 A surgical stainless steel. Brasseler USA's rigorous tempering and hardening process ensures outstanding sharpness and corrosion resistance.

▶ ANTERIOR SCALERS

25% OFF 6 OR MORE



SH6/H7	
#2, #4, #6	#8
\$32.25 ea.	\$34.25 ea.


Fine Point Offset Anterior Scaler
• 10° offset

SH5/48	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Fine Point Thin Hygienist Scaler/Hoe
• 2 mm hoe blade width

▶ POSTERIOR SCALERS

25% OFF 6 OR MORE



S204S	
#2, #4, #6	#8
\$32.25 ea.	\$34.25 ea.


Medium Sickle Scaler

S204SD	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Fine Posterior Sickle Scaler

▶ SAVANNAH SCALERS

25% OFF 6 OR MORE



SV1	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Savannah-1 - Anterior
• Shortened anterior sickle scaler end for interproximal calculus removal
• Disc end for removal of anterior lingual stain and calculus

SH5/SV1	
#6	
\$32.25 ea.	

Savannah-1 - Modified
• H5 anterior sickle scaler end for interproximal calculus removal
• Disc end for removal of anterior lingual stain and calculus

SV2	
#4, #6	
\$32.25 ea.	

Savannah-2 - Posterior
• Thin posterior sickle scaler designed for improved interproximal access
• Contra-angle design improves clinician comfort

SV3	
#6	#8
\$32.25 ea.	\$34.25 ea.

Savannah-3 - Posterior
• Thin contra-angle for better adaptation and enhanced posterior access

UNIVERSAL CURETTES

- Periodontal debridement instruments with rounded back, toe and two cutting edges.
- Can be used universally throughout the mouth.
- Manufactured from 440 A surgical stainless steel. Brasseler USA's rigorous tempering and hardening process ensures outstanding sharpness and corrosion resistance.

CURETTES

25% OFF 6 OR MORE



SBH5/6	
#2, #4, #6	#8
\$32.25 ea.	\$34.25 ea.

Barnhart Universal Curette
• For medium calculus removal



SC13/14	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Columbia Universal Curette with shortened blade and shank



SC4R/4L	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Columbia Universal Posterior Curette



SCR3/4	
#6	
\$32.25 ea.	

Ratcliff Universal Posterior Curette
• University of California



SL1/2	
#6	
\$32.25 ea.	

Langer Mandibular Posterior
• Same shank bend as SG11/12



SL3/4	
#6	
\$32.25 ea.	

Langer Maxillary Posterior
• Same shank bend as SG13/14



SL5/6	
#6	
\$32.25 ea.	

Langer Maxillary and Mandibular Anterior
• Same shank bend as SG5/6

SCURETTES

25% OFF 6 OR MORE

- Increase efficiency by adding S curettes to your tray setup.
- These unique instruments combine features of universal curettes and scalers in one instrument.



S137M	
#6	#8
\$32.25 ea.	\$34.25 ea.

Mini Long Universal Curette / Anterior Scaler



SN135	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Posterior Scaler with rounded back and pointed toe

GRACEY CURETTES

- Periodontal debridement instruments with a rounded back, toe and one cutting edge.
- Area-specific
- Multiple instruments are needed for instrumentation of the entire mouth.
- Manufactured from 440 A surgical stainless steel. Brasseler USA's rigorous tempering and hardening process ensures outstanding sharpness and corrosion resistance.

25% OFF 6 OR MORE



SG1/2	
#2, #4, #6	#8
\$32.25 ea.	\$34.25 ea.

Short contra-angle for incisors and cuspids



SG3/4	
#6	
\$32.25 ea.	

Short contra-angle for incisors and cuspids



SG5/6	
#6	#8
\$32.25 ea.	\$34.25 ea.

Medium contra-angle for anterior teeth and bicuspid



SG7/8	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

Medium contra-angle for buccal and lingual surfaces of bicuspid and molars



SG9/10	
#6	#8
\$32.25 ea.	\$34.25 ea.

Long contra-angle for buccal and lingual surfaces of molars and less accessible root surfaces



SG11/12	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

For mesial surfaces of bicuspid and molars



SG13/14	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

For distal surfaces of bicuspid and molars



SG15/16	
#6	#8
\$32.25 ea.	\$34.25 ea.

For improved access to mesial surfaces of posterior teeth
• Same shank bend as SG13/14



SG17/18	
#4, #6	#8
\$32.25 ea.	\$34.25 ea.

For improved access to distal surfaces of molars
• Shorter blade and longer terminal shank

All Brasseler Gracey Curettes are available in a variety of specialty patterns, such as Rigid, Long, Mini-Long and the new LTS™.



- Elongated terminal shank for deep access into pockets.
- 20% thinner blades than Mini-Long Graceys to reduce tissue distention and ease subgingival insertion.
- Blade length 50% shorter than Standard Graceys.
- Rigid shanks to reduce hand fatigue.
- Specially honed blade edges for razor sharpness.

IMPLANT MAINTENANCE INSTRUMENTS

▶ IMPLANTPRO™ TITANIUM IMPLANT CURETTES

25% OFF 6 OR MORE

- ImplantPro's 6Al-4V titanium is not heat treated, keeping the Rockwell C hardness at a low 25-31 HRC, and comparable to most implant abutments.
- Special blade finishing with a less-aggressive, polished edge ensures safe instrumentation around implants.
- Due to this unique blade finishing, hand sharpening is not recommended. Replacement tips are available.



TISL1/2

Titanium Langer 1/2

- Mandibular Posterior Universal Curette

\$89.95 ea.



TISL3/4

Titanium Langer 3/4

- Maxillary Posterior Universal Curette

\$89.95 ea.



TISL5/6

Titanium Langer 5/6

- Maxillary and Mandibular Anterior Universal Curette

\$89.95 ea.



TISC4R/4L

Titanium Columbia 4R/4L

- Posterior Universal Curette

\$89.95 ea.

ADVANCED HYGIENE INSTRUMENTS

▶ DIAMOND TIPPED FINISHING FILES

25% OFF 6 OR MORE

- 360° diamond coating in micro fine grit for maximum versatility and less-aggressive instrumentation



DTF1/2

Recommended for buccal and lingual furcations and root depressions

#6

- Long shank design improves access to deeper perio pockets with bifurcations and trifurcations

\$69.95 ea.




DTF3/4

Recommended for mesial root surfaces and distal depressions and line angles

#6

\$69.95 ea.



DTF3/4M

Recommended for mesial root surfaces and distal depressions and line angles

#6

- Modified F3 end perfectly complements anatomy
- Shortened shanks for enhanced clinician control

\$69.95 ea.

HYGIENEPRO™ ULTRASONIC SCALER INSERTS

24 Ergonomic Scaler Inserts to Provide Greater Comfort and Efficiency for You and Your Patients!



Wide variety of tip styles
available in both 25K and 30K frequency options to suit each user's needs

Competitive compatibility
with all magnetostrictive devices for ease and versatility

Ergonomic grip designs
feature textured non-slip surfaces for enhanced clinician comfort and reduced hand fatigue

Patented laminate solder joint
to ensure outstanding durability, consistent spray and uniform vibration



**#10
UNIVERSAL**
Available in slotted and direct focus water delivery options.



**#100
SLIM**
Available in slotted and direct focus water delivery options.



**#1000
TRIPLE BEND**
Available in slotted water delivery option only.



**#3
BEAVER TAIL**
Available in slotted water delivery option only.

▶	HARD GRIP	Slotted	\$139.99 ea.
		Direct Focus	\$157.99 ea.

▶	SOFT GRIP	Slotted	\$173.99 ea.
		Direct Focus	\$173.99 ea.



35% OFF 6 OR MORE!
AS LOW AS \$90.99 EACH!

Is it time to replace your ultrasonic scaler inserts?

Did you know that the last 1-2mm of an insert tip are the most active?
After 1mm of tip wear, 25% of the tip's efficiency is lost.
After 2mm of tip wear, 50% of the tip's efficiency is lost.

HYGIENEPRO™ AIR

Retail: \$599.00
Special: \$499.00

360° Swivel



5021318U0

Maximum Motor Speed: 2,600 rpm
Weight: 2.5 oz.
Length: 4.0 inches
Connects Directly to 4 Hole Tubing
1 Year Warranty

DISPOSABLE PROPHY ANGLES

BUY 2 BULK BOXES FOR \$999.00,
GET A HYGIENEPRO™ AIR MOTOR FREE



- Box of 100 or Bulk Box of 1200
- Latex Free

100 ct.		1200 ct.	
Firm	Soft	Firm	Soft
DPAF	DPAS	DPAF1200	DPAS1200
\$64.95/Box		\$697.25/Box	

Now Available
in Our New
Bulk Box!

DISPOSABLE PROPHY BRUSHES



- Box of 100
- Latex Free

DPAB
\$83.95/Box

DISPOSABLE PROPHY CUPS

30% OFF 3 OR MORE

SCREW TYPE



- Box of 144
- Latex Free

Firm	Soft
DPSCCF	DPSCCS
\$39.95/Box	

SNAP-ON TYPE



- Box of 144
- Latex Free

Firm	Soft
DPSNCF	DPSNCS
\$37.95/Box	

LATCH TYPE



- Box of 144
- Latex Free


Firm	Soft
DPLACF	DPLACS
\$48.95/Box	

GLISTEN™ PROPHY PASTE

30% OFF 3 OR MORE

Retail: \$42.75/Box

- Contains Fluoride
- Gluten Free
- Box contains 200 unit dose cups and 1 autoclavable prophy ring



		Fine	Medium	Coarse	Extra-Coarse	Kids Medium-Fine
Mint		PRP-MF	PRP-MM	PRP-MC	PRP-MXC	N/A
Cherry		N/A	PRP-CM	PRP-CC	N/A	N/A
Root Beer		N/A	PRP-RBM	PRP-RBC	N/A	N/A
Bubble Gum		N/A	PRP-BM	PRP-BC	N/A	N/A
Chocolate		N/A	PRP-CHOM	PRP-CHOC	N/A	N/A
Kids Asst. (Bubble Gum, Root Beer, Chocolate)		N/A	N/A	N/A	N/A	PRP01CH
Adult Asst. (Mint, Chocolate, Cherry)		N/A	PRP01AD	PRP02AD	N/A	N/A

FORZA™ V3 PIEZO ULTRASONIC UNIT

Elegant and Intuitive. These two words best describe the new Forza V3 piezo ultrasonic from Brasseler USA. At only 33 grams the slim, ergonomic design of the new Forza V3 handpiece provides excellent access while reducing hand fatigue. Bright LED illumination from dual fiber optics delivers clearer vision and maximum visibility. Equally important is the Piezo Accelerator mode that instinctively provides additional power at the moment you need it most. The compact, space saving unit complements the simple-to-operate interface and can be easily installed in any dental operator. With a wide range of piezo ultrasonic tips available, the Forza V3 is ideal for a variety of clinical procedures including general scaling, perio and endo.

- Dual LED Fiber Optics
- Slim & Balanced Handpiece Design
- Compact Control Unit
- Simple Operation Control Panel
- General, Perio & Endo Applications
- Easy Power & Water Adjustment

Actual Size



Featuring Piezo Accelerator Technology

5023152U0

Forza V3 LED Complete Set*

Retail: \$1,499.00

Special: \$1,199.00



5023153U0

Forza V3 Non-Optic Complete Set*

Retail: \$1,149.00

Special: \$919.00



* Set Contains: Control Unit with Handpiece Cord, Handpiece, Control Unit Holder, Handpiece Holder, Water Filter Set, Water Supply Connector, Foot Control, Tip Wrench with Torque Limiter, AC Adapter
Three Scaling Tips Included - G1, G9, G11 (\$282 Value)

► PIEZO ULTRASONIC TIPS 30% OFF 3 OR MORE

SCALING

G1  \$94.00

G6  \$94.00

G9  \$94.00

G11  \$94.00

PERIO

P1  \$142.00

P1D  \$152.00

P10  \$142.00

P20  \$142.00

P25R  \$142.00

P25L  \$142.00

P26R  \$142.00

P26L  \$142.00

BRASSELER USA

One Brasseler Blvd.
Savannah, Ga 31419

Toll Free: 800.841.4522

Phone: 912.925.8525

Fax: 888.610.1937

BrasselerUSA.com



BRASSELER
USA® BY YOUR SIDE

©2015 Brasseler USA. All rights reserved. Printed in the U.S.A.

B-4195-US-70M-04.15-PANA



BRASSELER
USA®

HYGIENE | 2015
PROMOTIONS